

Cell Growth Protocol for Jurkat Cell Line

From: HudsonAlpha/Caltech ENCODE group

Date: 8/27/08

Prepared by: Norma Neff and Tim Reddy

Jurkat (ATCC number TIB-152) cell culture and formaldehyde cross-linking

Jurkat Clone E6-1 is a human T lymphoblastoid cell line derived from an acute T cell leukemia. The cells are suspension lymphoblasts. The karyotype is pseudodiploid human male cell line. The modal chromosome number is 46, occurring in 74% with polyploidy at 5.3%.

Cell culture protocol:

Growth medium: Advanced RPMI 1640 (Gibco/Invitrogen) + 10% fetal bovine serum (Hyclone) + 10 mM Hepes + 100 units/ml penicillin + 100 µg/ml streptomycin + 5% CO₂ at 37°C.

Liquid Nitrogen Storage: Complete growth medium supplemented with 5% (v/v) DMSO in 1ml aliquots of approximately 5×10^6 cells.

1. Thaw a 1-ml aliquot of cells as quickly as possible in water bath at 37°C. Transfer cells to 9 ml warm media in a 15-ml conical tube. Mix gently. Centrifuge at 1,200 rpm for 5 minutes to pellet cells. Discard media and resuspend pellet gently in 10 ml warm medium. Divide cells into two T-25 flasks containing 5 ml warm media. Place in incubator. After two days, remove the media and add fresh media.
2. When cells are $6-8 \times 10^5$ cells/ml, split them 1:4 with fresh media. Add appropriate aliquots of the cell suspension to new culture vessels (T25 = 10 ml; T75 = 50 ml; T150 = 100 ml maximum volume). Grow cells to no more than 8×10^5 cells/ml. Disperse clumps gently for counting.

Large-scale cell growth (2-4 liters), cross-linking and harvest:

3. Grow cells in T150 flasks to a concentration of $5-8 \times 10^5$ cells/ml. Divide cells into four pools. Save aliquots of cells (1×10^6 to 5×10^7) from each pool for DNA or other types of analysis.
4. Fix cells as recommended for K562 cells. Concentrate cells by centrifugation at 1,200 rpm for 5 minutes at room temperature and resuspend in RPMI with no additives at 2×10^7 cells/ml. RPMI + 2% formaldehyde (freshly made) is added in equal volume to cells. Mix on rotator for 10 minutes at room temperature. Add glycine to 0.125 M, mix on a

rotator at room temperature for 5 minutes. Centrifuge at 800 rpm at room temperature for 10 minutes. Remove supernatant. Resuspend cells in cold PBS, pH 7.4. Transfer 10^8 cells to 15-ml conical tubes. Centrifuge at 800 rpm for 5 minutes at 4°C. Aspirate off PBS and freeze pellets on dry ice. Store at -80°C.